


# InLab<sup>®</sup> Redox Pt805

ORP half-cell with platinum ring

Scale 1:1



<b>Sensor type</b>	ORP half-cell
<b>Order number</b>	59904377
<b>Specifications</b>	
pH range	-
zero point (pH)	-
temperature range	-30...130 °C
temperature probe	-
shaft material	glass
membrane glass	-
membrane resistance	-
reference system	-
reference diaphragm	-
reference electrolyte	-
bridge electrolyte	-
electrode head	S7
cable length	-
connector	-
cell constant (cm <sup>-1</sup> )	-
cell type	-
indicating metal	platinum ring
<b>Storage</b>	
dry	